



LABORATOIRES PRODENE KLINT 8, rue Léon Jouhaux – ZA PARIEST CROISSY BEAUBOURG F - 77435 MARNE LA VALLEE CEDEX 2	SAFETY SHEET FINISHED PRODUCT	Approval: 
	File: C:/FICHE/FT	Edition: 0 / 12.02.09 Revision: 8 / 06.09.10 Page: 1/3

SHEET No.1200-FT



PREVEN`S PARIS LOTION

DEFINITION

PREVEN`S PARIS LOTION is an antiseptic* and hypoallergenic lotion for hygienic handwashing with rinsing.

PREVEN`S PARIS LOTION complies with the decree of 08/09/1999 and later laws concerning equipment cleaning products which are liable to come into contact with foodstuffs.

* formulated to minimise the risk of allergic reaction

COMPOSITION

Preparation made from non ionic surfactants, didecyltrimonium chloride CAS No.7173-51-5 (0.5% w/w), polyaminopropyl biguanide CAS No.27083-27-8 (0.4% w/w) and isopropyl alcohol CAS No.67-63-0 (0.2% w/w). Without betaine, cocamide dea, colorant and perfume


TESTS

1. *Organoleptic characteristics*

- *Aspect*: liquid clear slightly opalescent
- *Color*: colourless to pale yellow
- *Odour*: characteristic

2. *Physico-chemical characteristics (at 20°C)*

- *pH* (direct): 5.5 – 6.5
- *Density*: 1.004 – 1,0014
- *Viscosity* (Brookfield at 20°C): 1600 - 2600 mPa.s

LABORATOIRES PRODENE KLINT 8, rue Léon Jouhaux – ZA PARIEST CROISSY BEAUBOURG F - 77435 MARNE LA VALLEE CEDEX 2	SAFETY SHEET FINISHED PRODUCT	Approval: 
	File: C:/FICHE/FT	Edition: 0 / 12.02.09 Revision: 8 / 06.09.10 Page: 2/3

INSTRUCTIONS FOR USE:

Apply a 1,5 ml dose of product to wet skin, rub the hands, wrists, forearms, between the fingers and around the nails for at least 30 seconds; rinse carefully with water; dry using a single-use towel.

APPLICATION

- **Agri-food environment:** communal institutions, catering industry, production plant...
- **Medical environment:** private hospitals, hospitals, surgeries...

MICROBIOLOGICAL PROPERTIES

PREVEN`S PARIS LOTION complies with standards :

- EN 1040, 30s
- EN 1040, on *Listeria monocytogenes*, *Salmonella enterica* and *MRSA*
- NF T 72-170, state of dirtiness
- EN 1276, state of dirtiness
- EN 1499, 1,5mL, 30 seconds
- Pr EN 13727, clean conditions, 30s
- EN 1275, on *Candida albicans*
- NF T 72-301, *Candida albicans* in 30s
- According to the methodology of Standard EN 14476, 30s on *Adenovirus*, *Rotavirus*, *Poliovirus*, *HIV* and *H1N1*


PRECAUTIONS IN USE

Keep from freezing. Not to be used with other antiseptics or other products.

- **R 52/53** : Harmful to aquatic organisms, may cause long-term adverse effects in the aquatic environment.

STABILITY

In its original packaging and under normal storage conditions, **PREVEN`S PARIS LOTION** can be kept for 3 years.

LABORATOIRES PRODENE KLINT 8, rue Léon Jouhaux – ZA PARIEST CROISSY BEAUBOURG F - 77435 MARNE LA VALLEE CEDEX 2	SAFETY SHEET FINISHED PRODUCT	Approval: 
	<u>File:</u> C:/FICHE/FT	Edition: 0 / 12.02.09 Revision: 8 / 06.09.10 Page: 3/3

PACK SIZE

- 12 x 300mL with long-spouted pump
- 6 x 1L with long-spouted pump

Respect/Comfort for the user

Optimum skin tolerance with:
Formulation based on gentle surfactants
Dermatologically tested
Optimization of the dose and contact time
Restricted biocide active ingredients thanks to a synergy of action
Smooth to the touch

Respect for the environment

Limited environmental impact:
Free from EDTA (bioaccumulative), free from betaine (slow to biodegrade)
Contains natural and plant-based surfactants which are 100% biodegradable

Proven maximum effectiveness

Wide spectrum of action: bactericidal, levuricidal, virucidal
Tested in Cofrac (French Committee of Accreditation) accredited laboratories