



LABORATOIRES PRODENE KLINT 8, Rue Léon Jouhaux – Z.A. Pariest CROISSY BEAUBOURG F – 77435 MARNE LA VALLEE CEDEX 2	TECHNICAL INFORMATION FINISHED PRODUCT	Approval: 
	File : C:/FICHE/FT	Edition : 0 / 15.03.04 Révision : 12 / 21.09.2010 Page : 1/2

SHEET N°928-FT



DASR

DEFINITION

DASR is a very powerful ALDEHYDE-FREE and QUATERNARY AMMONIUM-FREE contact disinfectant, designed to disinfect by manual spraying, accredited by the Minister for Agriculture n°2030132, bactericidal treatment of 100 % animal and vegetable origin, no rinsing required, no need to wipe walls, surfaces and floors after application.

DASR : product complying with decree of 08/09/1999 and subsequent texts regarding equipment cleaning products possibly in contact with foodstuffs

COMPOSITION

Isopropyl alcohol CAS n°67-63-0 (30% m/m).

TESTS

1. Organoleptic characteristics

- **Appearance** : homogenous and limpid liquid
- **Odour** : alcohol
- **Colour** : colourless becoming yellow over time

2. Physico-chemical characteristics (at 20 °C)


- **pH (pure)**: 5,5 - 6,0

APPLICATIONS

For manual spraying. Spray **DASR** approximately 40 centimetres from the surface to be disinfected. Do not wipe.

Preparation and transport of feed for domestic animals.

Transportation, storage, industrial processing, sale of products of animal and vegetal origin.

LABORATOIRES PRODENE KLINT 8, Rue Léon Jouhaux – Z.A. Pariest CROISSY BEAUBOURG F – 77435 MARNE LA VALLEE CEDEX 2	TECHNICAL INFORMATION FINISHED PRODUCT	Approval: 
	File : C:/FICHE/FT	Edition : 0 / 15.03.04 Révision : 12 / 21.09.10 Page : 2/2

MICROBIOLOGICAL PROPERTIES:

DASR complies with the standards

- EN 1040, 5min
- EN 1276, clean conditions, 30s
- EN 1276, clean conditions, active on *Listeria monocytogenes* in 5min
- EN 1276, clean conditions active on salmonella
- NFT 72-281 bactericide, 15min
- NF EN 13697, 5min, *Pseudomonas aeruginosa*, *Escherichia coli*, *Staphylococcus aureus*, *Enterococcus hirae*
- NF EN 13697, 15min, *Candida albicans*
- EN 1650, *Candida albicans*, 1min
- EN 1275, *Candida albicans*, 15min
- According to the methodology of standard EN 14476+A1, active on H1N1 in 30s

PRECAUTIONS D'EMPLOI

Xi – IRRITANT

- **R 10** : Inflammable
- **R 36** : Irritant to the eyes.
- **R 67** : Vapours may cause drowsiness and dizziness.
- **S 23** : Do not breathe vapour
- **S 26** : In the event of contact with the eyes, wash immediately and copiously with water and consult a specialist.

PACKAGING

6 x 1L + 1 spray

6 x 750mL

2 x 5L

STORAGE

In its original sealed packaging, under normal storage conditions, **DASR** will keep for 3 years
May be damaged by freezing